

FocusDirect™ Secure Interface

*Seamless Connectivity with Your
Laboratory's Information System (LIS)*

FocusDirect™ is an interface that provides fully integrated electronic connectivity with Focus Diagnostics via your laboratory or hospital information system (LIS/HIS).

A FocusDirect™ interface enables a seamless two-way transmission of referral test orders and patient results using a secure Virtual Private Network (VPN) socket connection.

Benefits of Direct Interface

- Secure VPN socket connectivity providing a seamless secure tunnel protecting your patient's information
- Efficient transmission of original data eliminates transcription errors and double data entry
- User-definable communications to transmit at the laboratory's convenience
- Library of existing LIS protocols reduces time required for implementation
- Data transmission protocols are designed to industry standard HL7 (or as required)
- Technical support at each phase of implementation

Implementing FocusDirect™ Interface – Partner Team Approach

Once your laboratory is ready to implement a FocusDirect™ interface, our goal is to streamline the process using a team approach. Focus will assess your laboratory system requirements and dedicate a team to implement a direct two-way (orders and results) interface between both locations. Your Focus Project Manager will oversee every detail of the interface plan, working as a team leader directly with the client and information system vendor to facilitate the design, testing and final implementation.

FocusDirect™ Interfaces are currently implemented on numerous laboratory information systems. Contact your Focus Sales Account Manager for a complete list of current LIS vendor interfaces.

All Focus connectivity solutions are designed to comply with HIPAA regulations.

Please contact your Focus Account Manager to learn more about FocusLink® or other connectivity options from Focus Diagnostics.

800-445-0185 or **714-220-1900**

www.focusdx.com

Focus Diagnostics
Cypress, California USA

©Copyright Focus Diagnostics, Inc. 2009
RLCFD0309



*Web-based or dial-up
connection directly to
Focus LIS allows quick
look-up anytime.*

*FocusConnect™ Suite of
Connectivity Solutions*

The source for infectious disease
diagnostic solutions

